



上海源叶生物科技有限公司
Shanghai yuanye Bio-Technology Co., Ltd
电话: 021-61312973 传真: 021-55068248
网址: www.shyuanye.com
邮箱: shyysw@sina.com

谷胱甘肽还原酶酶活测试

PRINCIPLE:



Abbreviations used:

β -NADPH = β -Nicotinamide Adenine Dinucleotide Phosphate, Reduced Form

β -NADP = β -Nicotinamide Adenine Dinucleotide Phosphate, Oxidized Form

GSSG = Glutathione, Oxidized Form

GSH = Glutathione, Reduced Form

CONDITIONS: T = 25°C, pH = 7.6, A_{340nm}, Light Path = 1 cm

METHOD: Continuous Spectrophotometric Rate Determination

REAGENTS:

A. 100 mM Potassium Phosphate Buffer with 3.4 mM Ethylenediaminetetraacetic

Acid (EDTA), pH 7.6 at 25°C (Prepare 200 ml in deionized water using Potassium

Phosphate, Monobasic, Anhydrous, and Ethylenediaminetetraacetic Acid,

Dipotassium Salt, Adjust to pH 7.6 at 25°C with 1 M KOH.)

B. 30 mM Glutathione Substrate Solution (GSSG)

(Prepare 5 ml in deionized water using Glutathione, Oxidized Form, Disodium Salt.)

C. 0.8 mM β -Nicotinamide Adenine Dinucleotide Phosphate, Reduced Form

Solution (β -NADPH) (Prepare 5 ml in cold Reagent A using β -Nicotinamide Adenine

Dinucleotide Phosphate, Tetrasodium Salt.)

D. 1.0% (w/v) Bovine Serum Albumin (BSA)

(Prepare 100 ml in Reagent A using Albumin, Bovine. This solution should be kept cold.)

E. Glutathione Reductase Enzyme Solution



上海源叶生物科技有限公司
Shanghai yuanye Bio-Technology Co., Ltd
电话: 021-61312973 传真: 021-55068248
网址: www.shyuanye.com
邮箱: shyysw@sina.com

(Immediately before use, prepare a solution containing 0.30 - 0.60 unit/ml of

Glutathione Reductase in cold Reagent D.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes:

	Test	Blank
Deionized Water	0.65	0.65
Reagent A (Buffer)	1.50	1.50
Reagent B (GSSG)	0.10	0.10
Reagent C (β -NADPH)	0.35	0.35
Reagent D (BSA)	0.30	0.40

Mix by inversion and equilibrate to 25°C. Monitor the $A_{340\text{nm}}$ until constant, using a suitably thermostatted spectrophotometer. Then add:

Reagent E (Enzyme Solution)	0.10	-----
-----------------------------	------	-------

Immediately mix by inversion and record the decrease in the $A_{340\text{nm}}$ for approximately 5 minutes. Obtain the $\Delta A_{340\text{nm}}/\text{minute}$ using the maximum linear rate for both the Test and Blank.

CALCULATIONS:

$$\text{Units/ml enzyme} = \frac{(\Delta A_{340\text{nm}}/\text{min Test} - \Delta A_{340\text{nm}}/\text{min Blank})(3)(\text{df})}{(6.22)(0.1)}$$

3 = Total volume (in milliliters) of assay

df = Dilution factor

6.22 = Millimolar extinction coefficient of β -NADPH at 340 nm

0.1 = Volume (in milliliters) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

UNIT DEFINITION:

One unit will reduce 1.0 μmole of oxidized glutathione per minute at pH 7.6 at 25°C.

FINAL ASSAY CONCENTRATION:

In a 3.00 ml reaction mix, the final concentrations are 75 mM potassium phosphate, 2.6 mM ethylenediaminetetraacetic acid, 1 mM glutathione, 0.09 mM β -nicotinamide adenine dinucleotide phosphate, reduced form, 0.13% (w/v) bovine serum albumin, and 0.03 - 0.06 unit of glutathione reductase.