


MASS SPECTROMETRY REPORT



Sample Description		Instrument		SHIMADZU LCMS-2010	
Analyzed date: 2015-01-08		Probe:	ESI	Probe Bias:	+4.5kv
Analyst:	QIU	Nebulizer Gas Flow:	1.5L/min	Detector:	1.5kv
Sample:	Angiotensin II CAS#4474-91-3	CDL:	-20.0v	T. Flow:	0.2ml/min
M.W.:	1046.2	CDL Temp.:	250 °C	B. Conc.:	50% H2O/50% ACN
Lot. No.:	P150112-LR193974	Block Temp.:	200 °C		