



上海源叶生物科技有限公司
Shanghai yuanye Bio-Technology Co., Ltd
电话: 021-61312973 传真: 021-55068248
网址: www.shyuanye.com
邮箱: shyysw@sina.com

丁酰胆碱酯酶活性测试

PRINCIPLE:

Butyrylcholine + H₂O Cholinesterase, Butyryl> Butyrate + Choline

CONDITIONS: T = 37° C, pH = 8.0

METHOD: Titrimetric

REAGENTS:

A. 1600 mM Magnesium Chloride Solution

(Prepare 15 ml in deionized water using Magnesium Chloride, Hexahydrate,.)

B. 1000 mM Sodium Chloride Solution

(Prepare 50 ml in deionized water using Sodium Chloride.)

C. 4 mM Butyrylcholine Chloride Solution (BuCholine)

(Prepare by dissolving 336 mg of Butyrylcholine Chloride,, in 350 ml deionized water. Then add 10 ml Reagent A, and 40 ml Reagent B for a total volume of 400 ml.)

D. 20 mM Sodium Hydroxide Solution-Standardized (NaOH)

(Prepare 50 ml in cold deionized water using Sodium Hydroxide,)

E. Cholinesterase, Butyryl Enzyme Solution

(Immediately before use, prepare a solution containing 30 - 60 units/ml of Cholinesterase, Butyryl in cold deionized water.)

PROCEDURE:

Using a suitable pH meter in conjunction with a magnetic stirrer, pipette (in milliliters) the following

reagents into a suitably thermostatted titration vessel:

Test

Reagent C (BuCholine) 50.00

Equilibrate to 37° C. Adjust to pH 8.5 with Reagent D (NaOH), using a burette. Then add:

Reagent E (Enzyme Solution) 0.40

Run the reaction for 1 - 5 minutes. Record the time when the pH reaches 8.0. Maintain the pH of the reaction mix at pH 8.0 by the addition of small volumes (0.05 ml) of Reagent D. Record the volume of Reagent D used to maintain the pH at 8.0 and the time required

CALCULATION:



上海源叶生物科技有限公司
Shanghai yuanye Bio-Technology Co., Ltd
电话: 021-61312973 传真: 021-55068248
网址: www.shyuanye.com
邮箱: shyysw@sina.com

$$\text{Units/mg enzyme} = \frac{(\text{Molarity of NaOH}) (\text{NaOH}) (1000)}{(T) (\text{mg enzyme/RM})}$$

NaOH = Volume (in milliliters) of Reagent D used in the assay

1000 = Conversion from millimoles to micromoles (Unit definition)

T = Time of assay

RM = Reaction Mix

UNIT DEFINITION:

One unit will hydrolyze 1.0 PROHRIbutyrylcholine to choline and butyrate per minute at pH 8.0 at

37° C.

INITIAL ASSAY CONCENTRATIONS:

In a 50.4 ml reaction mix, the initial concentrations are 40 mM magnesium chloride, 99 mM sodium chloride, 4 mM butyrylcholine chloride and 12 - 24 units cholinesterase, butyryl.