

MASS SPECTROMETRY REPORT



Sample Description		Instrument		SHIMADZU LCMS-2020	
Analyzed date:	2018-9-26	Probe:	ESI	Probe Bias:	+4.5kv
Analyst:	SHAO	Nebulizer Gas Flow:	1.5L/min	Detector:	1.5kv
Sample:	黑素细胞刺激素	CDL:	-20.0v	T. Flow:	0.2ml/min
M.W.:	1664.92	CDL Temp.:	250℃	B. Conc.:	50%H2O/50%ACN
Lot. No.:		Block Temp.:	200℃		