



上海源叶生物科技有限公司  
Shanghai yuanye Bio-Technology Co., Ltd  
电话: 021-61312973 传真: 021-55068248  
网址: www.shyuanye.com  
邮箱: shyysw@sina.com

产品名称: **Nortadalafil**  
产品别名: 去甲基他达那非; **Demethyl Tadalafil**

生物活性:					
<b>Description</b>	Nortadalafil is demethyl Tadalafil, which is a PDE5 inhibitor, currently marketed in pill form for treating erectile dysfunction (ED) under the name Cialis; and under the name Adcirca for the treatment of pulmonary arterial hypertension. IC50 value: Target:				
<b>Solvent&amp;Solubility</b>	<b>In Vitro:</b> <b>DMSO : 50 mg/mL (133.20 mM; Need ultrasonic)</b> <b>H<sub>2</sub>O : &lt; 0.1 mg/mL (insoluble)</b>				
		<b>Solvent Mass</b>	<b>1 mg</b>	<b>5 mg</b>	<b>10 mg</b>
		<b>Concentration</b>			
	<b>Preparing</b>	1 mM	2.6640 mL	13.3198 mL	26.6397 mL
	<b>Stock Solutions</b>	5 mM	0.5328 mL	2.6640 mL	5.3279 mL
		10 mM	0.2664 mL	1.3320 mL	2.6640 mL
<p>*请根据产品在不同溶剂中的溶解度选择合适的溶剂配制储备液。一旦配成溶液, 请分装保存, 避免反复冻融造成的产品失效。</p> <p>储备液的保存方式和期限: -80°C, 6 months; -20°C, 1 month。-80°C 储存时, 请在 6 个月内使用, -20°C 储存时, 请在 1 个月内使用。</p> <p><b>In Vivo:</b></p> <p>请根据您的实验动物和给药方式选择适当的溶解方案。以下溶解方案都请先按照 In Vitro 方式配制澄清的储备液, 再依次添加助溶剂:</p> <p>——为保证实验结果的可靠性, 澄清的储备液可以根据储存条件, 适当保存; 体内实验的工作液, 建议您现用现配, 当天使用; 以下溶剂前显示的百分比是指该溶剂在您配制终溶液中的体积占比; 如在配制过程中出现沉淀、析出现象, 可以通过加热和/或超声的方式助溶</p> <p>1.请依序添加每种溶剂: 10% DMSO→40% PEG300 →5% Tween-80 → 45% saline Solubility: ≥ 2.5 mg/mL (6.66 mM); Clear solution 此方案可获得 ≥ 2.5 mg/mL (6.66 mM, 饱和度未知) 的澄清溶液。 以 1 mL 工作液为例, 取 100 μL 25.0 mg/mL 的澄清 DMSO 储备液加到 400 μL PEG300 中, 混合均匀 向上述体系中加入 50 μL Tween-80, 混合均匀; 然后继续加入 450 μL 生理盐水定容至 1 mL。</p> <p>2.请依序添加每种溶剂: 10% DMSO →90% corn oil Solubility: ≥ 2.5 mg/mL (6.66 mM); Clear solution 此方案可获得 ≥ 2.5 mg/mL (6.66 mM, 饱和度未知) 的澄清溶液, 此方案不适用于实验周期在半个月以上的实验。 以 1 mL 工作液为例, 取 100 μL 25.0 mg/mL 的澄清 DMSO 储备液加到 900 μL 玉米油中, 混合均匀。</p>					
<p>[1]. Beghyn, et al. PDE5 inhibitors: An original access to novel potent arylated analogues of tadalafil. Bioorganic &amp; Medicinal Chemistry Letters (2007), 17(3), 789-792.</p> <p>[2]. Daugan, et al. The discovery of tadalafil: a novel and highly selective PDE5 inhibitor. 2:</p>					



上海源叶生物科技有限公司  
Shanghai yuanye Bio-Technology Co., Ltd  
电话: 021-61312973 传真: 021-55068248  
网址: [www.shyuanye.com](http://www.shyuanye.com)  
邮箱: [shyysw@sina.com](mailto:shyysw@sina.com)

#### References

- 2,3,6,7,12,12a-hexahydropyrazino[1',2':1,6]pyrido[3,4-b]indole-1,4-dione analogue. Journal of Medicinal Chemistry (2003), 46(21), 4533-4542.
- [3]. Porst H, et al. Chronic PDE-5 inhibition in patients with erectile dysfunction - a treatment approach using tadalafil once-daily. Expert Opin Pharmacother. 2012 Jul;13(10):1481-94.
- [4]. Frey MK, Lang I. Tadalafil for the treatment of pulmonary arterial hypertension. Expert Opin Pharmacother. 2012 Apr;13(5):747-55.



# 源叶生物